



# Leopard Spinal Brace





## Dual Coding

Most BDO spinal braces are dual-coded by PDAC. That means that are designed to be dispensed as is (OTS) or to be custom fitted by an experienced professional.

The polymer plates may be heat remolded to conform to the torso. For other custom modifications, see [www.bluedortho.com](http://www.bluedortho.com).

## Satisfaction

When you fit a Blue Diamond Orthopedic brace, you'll be rewarded with the gratitude of trusting patients and the respect of referring physicians.

## Ask Us About

- Custom-branding
- Becoming a sales representative

## Warranty

Blue Diamond Orthopedic warrants its spinal braces for one year.

# LEOPARD SP

## AP Control LSO

The Universal Leopard SP spinal brace is indicated for the patient with moderate to severe low-back pain. The luxurious tri-laminate fabric and the rigid Kydex posterior plate combine to provide comfortable lumbo-sacral control.

The quad-strap compression system snugs the brace firmly.

The Leopard SP is easy to fit. Set the size by trimming and reattaching the lateral extensions, reset the straps, don and tighten. That's it!

PDAC Approved: L0631/L0648



The Best Value in Spinal Bracing

800-867-6070 | [www.bluedortho.com](http://www.bluedortho.com)

# LEOPARD SP



## SPECIFICATIONS

- AP Control LSO
- Plush tri-laminate fabric
- Superior sagittal control
- One size – easy to stock
- Conforms to any body shape
- Adjustable quad-pull compression

## UNIVERSAL SIZE

One size fits 30" to 65"

## PDAC APPROVED

- L0631 Custom Fit
- L0648 Off the Shelf (OTS)

## PRECISION ADJUSTMENT

The Leopard SP's genuine Kydex anterior and posterior control plates can be reformed with a heat gun for a precise fit.



Blue Diamond Orthopedic | 6439 Milner Blvd, Suite 4, Orlando, FL 32809 | 800-867-6070 | 407-613-2001 | [bluedortho.com](http://bluedortho.com)

© 2019 Blue Diamond Orthopedic. All rights reserved. Patent Pending

20200625

# LEOPARD SL

## APL Control LSO

The Universal Leopard SL spinal brace is indicated for the patient with moderate to severe low-back pain. The breathable mesh and the rigid Kydex posterior and lateral plates combine to provide bi-planar LSO control.

The quad-strap compression system snugs the brace firmly.

The Leopard SL is easy to fit. Set the size by trimming and reattaching the lateral extensions, reset the straps, affix the lateral control pouches, don and tighten. That's it!

PDAC Approved: L0637/L0650



The Best Value in Spinal Bracing

800-867-6070 | [www.bluedortho.com](http://www.bluedortho.com)



# LEOPARD SL



## SPECIFICATIONS

- APL Control LSO
- Plush tri-laminate fabric
- Superior sagittal control
- One size – easy to stock
- Conforms to any body shape
- Lateral plates for coronal control
- Adjustable quad-pull compression

## UNIVERSAL SIZE

One size fits 30" to 65"

## PDAC APPROVED

- L0637 Custom Fit
- L0650 Off the Shelf (OTS)

## PRECISION ADJUSTMENT

The Leopard SL's genuine Kydex anterior, posterior, and lateral control plates can be reformed with a heat gun for a precise fit.



Blue Diamond Orthopedic | 6439 Milner Blvd, Suite 4, Orlando, FL 32809 | 800-867-6070 | 407-613-2001 | [bluedortho.com](http://bluedortho.com)

© 2019 Blue Diamond Orthopedic. All rights reserved. Patent Pending

20200625

# LEOPARD TL

## Extension Assist TLSO

The Leopard TL spinal brace is an excellent TLSO for the treatment of compression fractures, excessive kyphosis, or for post-fusion stabilization.

The quad-strap compression system snugs the brace firmly.

The Leopard TL is easy to fit. Set the size by trimming and reattaching the lateral extensions, reset the straps, adjust the shoulder straps, don and tighten.

That's it!

PDAC Approved: L0456/L0457



The Best Value in Spinal Bracing

800-867-6070 | [www.bluedortho.com](http://www.bluedortho.com)

# LEOPARD TL



## SPECIFICATIONS

- Extension Assist TLSO
- Plush tri-laminate fabric
- One size – easy to stock
- Conforms to any body shape
- Adjustable-height posterior strut
- Adjustable shoulder/axilla straps
- Sacral to thoracic flexion control
- Adjustable quad-pull compression

## PDAC APPROVED

- L0456 Custom Fit
- L0457 Off the Shelf (OTS)

## PRECISION ADJUSTMENT

The Leopard TL's genuine Kydex anterior and posterior control plates can be reformed with a heat gun for a precise fit.

## UNIVERSAL SIZE

One size fits 30" to 65"



Blue Diamond Orthopedic | 6439 Milner Blvd, Suite 4, Orlando, FL 32809 | 800-867-6070 | 407-613-2001 | [bluedortho.com](http://bluedortho.com)

© 2019 Blue Diamond Orthopedic. All rights reserved. Patent Pending

20200625



# LEOPARD TLSSL

## Quad Control TLSO



The Leopard TLSSL spinal brace delivers maximum spinal control. The posterior, anterior and long bi-lateral control plates overlap to provide circumferential support & control.

With its sternal STEP and attached shoulder straps, the Leopard TLSSL offers triplanar control from the symphysis pubis to the sternal notch anteriorly and from the gluteal fold to the scapular spines posteriorly.

The quad-strap compression system snugs the brace firmly.

To fit the Leopard TLSSL, set the torso circumference by trimming and reattaching the lateral extensions, then replace the straps. Adjust the long lateral pouch and the shoulder straps, then don the brace, snug it to the torso, and set the STEP height.

PDAC Approved: L0464



The Best Value in Spinal Bracing

800-867-6070 | [www.bluedortho.com](http://www.bluedortho.com)

# LEOPARD TLSL



## SPECIFICATIONS

- Quad Control TLSO
- Posterior-lateral shell
- Plush tri-laminate fabric
- Smooth quad-pull compression
- Sacral to thoracic spinal control
- Adjustable shoulder/axilla straps
- Choice of STEP
  - I – Greater vertical range
  - II – Quick-adjust extension

## PDAC APPROVED

- L0464 Custom Fit

## PRECISION ADJUSTMENT

The Leopard TLSL anterior plate, posterior-lateral shell, lateral control plates, and sternal STEP plate can be reformed with a heat gun for a precise fit.

## UNIVERSAL SIZE

One size fits 30" to 65"



Blue Diamond Orthopedic | 6439 Milner Blvd, Suite 4, Orlando, FL 32809 | 800-867-6070 | 407-613-2001 | [bluedortho.com](http://bluedortho.com)

© 2019 Blue Diamond Orthopedic. All rights reserved. Patent Pending

20200625