

Panther MP 627



Moderate Control LSO

The Panther MP spinal brace is indicated for the patient with moderate low-back pain. The breathable mesh and rigid Kydex posterior plate combine to provide mild LSO control.

The string-and-pulley compression system allows the brace to be easily snugged, even by elderly arthritic hands.

The Panther MP is universal. It requires no tools to size. There's no need to cut or fold a belt. Just add or remove belt segments.

PDAC Approved: L0627/L0642



SCAN FOR
PANTHER CATALOG



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

Panther MP 627



SPECIFICATIONS

- Bilateral pull handles
- Moderate AP control
- Breathable-mesh belt
- Low-profile posterior plate
- Conforms to any body shape
- Smooth dual-pull compression

UNIVERSAL SIZE

- One size fits all.
- 25" to 66" and beyond! Just add or remove belt segments.

PDAC APPROVED

- L0627 Custom Fit
- L0642 Off the Shelf (OTS)

LORDOSIS

The posterior polymer plate is preset at 15°. It may be heat molded for exact fit during custom fitting.

