

Panther SP 6-31



AP Control LSO

The Panther SP spinal brace is indicated for the patient with moderate to severe low-back pain. The breathable mesh and the rigid Kydex posterior plate combine to provide comfortable LSO control.

The string-and-pulley compression system allows the brace to be easily snugged, even by elderly arthritic hands.

The Panther SP is universal. It requires no tools to size. There's no need to cut or fold a belt. Just add or remove belt segments.

PDAC Approved: L0631/L0648



SCAN FOR
PANTHER CATALOG



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

Panther SP 631



SPECIFICATIONS

- Bilateral pull handles
- Breathable-mesh belt
- Superior sagittal control
- Full-height posterior plate
- Conforms to any body shape
- Smooth dual-pull compression

UNIVERSAL SIZE

- One size fits all.
- 25" to 66" and beyond! Just add or remove belt segments.

PDAC APPROVED

- L0631 Custom Fit
- L0648 Off the Shelf (OTS)

LORDOSIS

The posterior polymer plate is preset at 15°. It may be heat molded for exact fit during custom fitting.

