



2025 WHOLESALE CATALOG

TABLE OF CONTENTS

- About BioSkin
- Knee
- Spine
- Foot & Ankle
- Upper Extremity
- Lower Extremity





pro.bioskin.com

CALL: 800-541-2455 | FAX: 541-482-2119 | EMAIL: orders@bioskin.com



INTERNATIONAL

We have distributors located all over the world.



A PARTNERSHIP FOR EXCEPTIONAL CARE

We know your goal is to provide exceptional care to your patients. Here are some of the ways our products help you do that:

- Thoughtfully-designed products prevent injury, reduce pain, and restore mobility
- Hypoallergenic materials are dry, cool, and comfortable for all-day use
- A wide range of sizes and design options provide the best fit possible
- Quality components are lightweight and durable for comfortable products that last
- Attention to details like carefully placed seams and non-constricting straps increase patient compliance






BIOSKIN™

A close-up, low-angle shot of a person's lower body and legs as they run on a green athletic track. The person is wearing a red shirt, black shorts, and a black knee brace on their right knee. Their right leg is extended forward, and their left leg is pushing off the ground. They are wearing blue and yellow running shoes. The background is a blurred green field and a chain-link fence under a clear blue sky.

KNEE BRACES

Product Shown: Crossfire OA

KNEE SIZE CHART

Contact Us to Place an Order

- Call: 800-541-2455
- Email: orders@bioskin.com
- Online: www.pro.bioskin.com

KNEE PRODUCTS					
Circumference with knee fully extended					
Size	*	Inches		Centimeters	
		7" above mid-patella	6" below mid-patella	18 cm above mid-patella	15 cm below mid-patella
XS	0	16 – 18	12 – 13.5	41 – 46	30 – 34
S	1	18 – 20	13.5 – 15	46 – 51	34 – 38
M	2	20 – 22	15 – 16.5	51 – 56	38 – 42
L	3	22 – 25	16.5 – 18.5	56 – 64	42 – 47
XL	4	25 – 28	18.5 – 21	64 – 71	47 – 53
XXL	5	28 – 32	21 – 24	71 – 81	53 – 61
3XL	6	32 – 36	24 – 27	81 – 91	61 – 69
4XL	7	36 – 40	27 – 30	91 – 102	69 – 76
5XL	8	40 – 44	30 – 33	102 – 112	76 – 84

Q BABY			
Circumference under patella with knee fully extended			
Size	*	Inches	Centimeters
XS - M	2	12 – 16	31 – 41
L - XXL	5	16 – 21	41 – 54

Please contact us for any custom sizing needs.

Product#	Size#	To indicate size, replace * with corresponding number from the sizing chart.
XXXX	*	

EVERY DETAIL MATTERS

A brace can't help if the patient won't wear it, so we obsess over details to make sure our products are both supportive and comfortable.

PROPRIETARY MATERIALS

- Medical-grade compression reduces pain and swelling
- Moisture activates the SkinLok™ feature and prevents migration

BUILT FOR COMFORT

- Elastic straps allow muscle expansion and minimize migration
- Die-cut velcro tabs eliminate sharp corners
- No seams behind the knee

QUALITY HINGES

- High quality hinges are both durable and lightweight



Lightweight



Dynamic Stretch



Breathable



Hypoallergenic

PERFECT FIT

Every leg is different and one-size-fits-all braces sacrifice efficacy and compliance. That's why we have multiple configurations as well as extended sizing from XS - 5XL.



WRAPAROUND DESIGN

Many of our braces are available with our unique wraparound design. The entire front panel is Velcro compatible, making it easy to achieve a perfect fit on virtually any leg. This is also a great option for patients with limited flexibility or hand strength.

PULL-ON DESIGN

Our pull-on knee braces are made with a highly compressive lycra-based fabric to maximize pain relief, control swelling, and enhance proprioception.



CROSSFIRE™

Any brace with CrossFire in the name incorporates innovative design features to maximize the therapeutic benefits of compression and increase comfort.



CROSSFIRE HYBRID WRAPAROUND DESIGN

An inner sleeve facilitates donning by holding the brace in place so patients have both hands available to adjust the fit. The sleeve provides compression to the entire knee to control swelling and deliver pain relief while the wraparound component of the brace also provides more sizing flexibility and a perfect fit.



ULTRA LIGHTWEIGHT MATERIAL BEHIND THE KNEE

CrossFire wraparound and pull-on braces feature our exceptionally lightweight Cirrus material behind the knee. This back panel is sewn in a slightly flexed position, virtually eliminating bunching behind the knee.



CROSSFIRE PULL-ON DESIGN

Some CrossFire braces are also available in a pull-on version that features the same lightweight and compressive material over the patella as shown in the wraparound version's donning sleeve. These braces also feature a more generous cut above the knee to accommodate a wider variety of leg shapes.

CROSSFIRE OA

DYNAMIC UNLOADER

L1851/L1843

Offloads the medial or lateral compartment to reduce osteoarthritis or meniscus pain. Also addresses patellofemoral pain.

- Spiral strap works in tandem with the hinge, calf, and thigh to create a three-point leverage system that opens the joint during extension. The strap relaxes during flexion
- Magnetic buckle attachment on the spiral strap allows for quick and easy application and removal of the brace
- Optional calf and thigh straps provide additional support
- Hybrid version features a pull on calf and wraparound thigh closure

PRODUCT	PART #
CrossFire OA	5190* (L or R)

Note: For a left/medial, right/lateral brace, order L.
For a right/medial, left/lateral, order R.



Optional patella traction strap (included) can adjust the tension on the patella to correct maltracking and relieve pain.

ROM HINGES: KNEE STABILIZATION WITH MOTION CONTROL

These braces provide high-level stabilization and limit flexion and extension with adjustable range-of-motion hinges.



GLADIATOR™ L1833/L1832

Our most supportive hinged knee brace.

- Non-elastic straps for maximum stabilization
- Rigid, aluminum, 15" uprights may be bent for a customized fit (L1832)
- Gladiator hinge: ROM hinge with flexion/extension stops at 15°, 30°, 60°, and 90°

PRODUCT	PART #
Gladiator Wraparound	5136*



CROSSFIRE ROM L1833/L1832

Exceptionally comfortable with a lighter weight hinge.

- Rigid, aluminum uprights are flared to contour to the leg's natural shape and may be bent for additional customization (L1832)
- SlimLok hinge: Lightweight, adjustable ROM hinge with three flexion and extension stops at 0°, 15°, and 30°
- Pull-on version includes a silicone gel pad around the patella to reduce swelling and decrease knee pain

PRODUCT	PART #
CrossFire ROM Wraparound	5180*
CrossFire ROM Pull-on with Gel	6180*



Q LOK™ ROM L1833/L1832

Patellar control in a highly supportive hinged knee brace.

- Traction strap controls patellar glide and tilt. Fits left or right leg
- Non-elastic straps control medial/lateral and anterior/posterior movement of the joint
- Rigid, aluminum uprights are flared to contour to the leg's natural shape and can be bent for additional customization (L1832)
- SlimLok hinge: Lightweight, adjustable ROM hinge with three flexion and extension stops at 0°, 15°, and 30°

PRODUCT	SPECIAL ORDER, LEAD TIMES VARY	PART #
Q Lok ROM Wraparound (interior donning sleeve)		5166*
Q Lok ROM Pull-on		6166*

LIGHTWEIGHT ALUMINUM HINGES: MODERATE STABILIZATION

Bicentric, aluminum hinges are lightweight yet supportive. Made of aircraft-grade aluminum for its unique strength-to-weight ratio.



CROSSFIRE ALUMINUM L1820

Made of aircraft-grade aluminum for its high strength to weight ration, these bicentric, aluminum hinges are lightweight and supportive.

- Inner donning sleeve for ease of application and light compression over the patella
- Exceptionally lightweight material is sewn in tension behind the knee to prevent uncomfortable bunching

PRODUCT

PART

CrossFire Aluminum | Wraparound

5186*



Q LOK™ DYNAMIC PATELLA TRACTION BRACE L1820

Heavy-duty pain relief and patella traction control from a lightweight knee brace.

- Closed patella design and elastic straps help to seat the patella in the trochlear groove
- Dynamic traction strap applies pressure on the patella as the knee bends

PRODUCT

PART

Q Lok Dynamic Patella Traction Brace - Stratus

6164*



HINGED KNEE BRACE L1820

Lightweight, hinged knee brace for general joint stabilization.

- Bicentric, aluminum hinges support ligaments and allow natural motion of the knee
- Pull-on version includes foam buttress around the patella for additional stabilization

PRODUCT

PART

Hinged Knee Brace | Wraparound

5133*

Hinged Knee Brace | Pull-on

4130*

CONFORMA™ HINGE: FLEXIBLE, LIGHTWEIGHT SUPPORT

The Conforma Hinge is made of lightweight, flexible nylon that conforms to the shape of any leg. Surprisingly lightweight, comfortable, and durable.



Q BRACE CROSSFIRE WITH CONFORMA HINGE L1820

A lightweight hinged knee brace with optional patella tracking strap.

- Lightweight, flexible hinges provide sagittal control without sacrificing comfort
- Universal design (Q strap works on either knee and can track in multiple directions)
- Prevents subluxation when worn with patella tracking strap
- Also available without a hinge

PRODUCT	PART #
Q Brace CrossFire with Conforma	5138*
Q Brace Pull-On (without hinge)	4161*
Q Brace Wraparound (no hinge)	5163*



CONFORMA™ HINGE: FLEXIBLE, LIGHTWEIGHT SUPPORT

The Conforma hinge is made of lightweight, flexible nylon that conforms to the shape of any leg. Surprisingly lightweight, comfortable, and durable.



VISCO CROSSFIRE WITH CONFORMA HINGE L1820

An all-purpose knee brace with a silicone gel ring perfect for relieving general knee pain and swelling.

- Soft, silicone gel ring around the patella reduces swelling and increases compression around the joint
- Also available without a hinge

PRODUCT	PART #
Visco CrossFire with Conforma	6137*
Visco Knee Sleeve with Straps (No hinge)	6171*



UNIVERSAL PREMIUM J™ WITH CONFORMA HINGE L1820

A comfortable solution for patellofemoral pain and joint stabilization.

- J-shaped silicon buttress forms a lateral wall to prevent subluxation and encourage proper alignment. Packaged with bilateral buttresses; simply remove one side to fit the left or right leg
- Elastic, anterior straps control the patella to correct glide and tilt. May be positioned on the right or left
- Also available left or right specific (6138* L or R), and without a hinge

PRODUCT	PART #
Universal Premium J w/ Conforma Hinge (Special order, lead times vary)	6138*U
Premium J with Conforma Hinge	6138* (L or R)



PATELLA STABILIZER WITH CONFORMA HINGE L1820

Stabilizes the patella, reduces swelling, and relieves anterior knee pain.

- A foam buttress encircles the patella providing passive massage, edema control, and stabilization

PRODUCT	PART #
Patella Stabilizer with Conforma Hinge	4132*

COMPRESSION PRODUCTS

Medical-grade compression provides pain relief.



PATELLA STABILIZER SLEEVE

Stabilizes the patella, reduces swelling, and relieves anterior knee pain.

- A foam buttress encircles the patella providing passive massage, edema control, and stabilization

PRODUCT	PART #
Patella Stabilizer Sleeve	4125*



ULTIMA KNEE SLEEVE™

A compressive knee sleeve for comfortable support and pain relief.

- Flexible stays on each side provide light stabilization to the knee
- Compression around the knee reduces pain and swelling

PRODUCT	PART #
Ultima Knee Sleeve	4110*O



Q BABY™

An effective patellar tendon strap with a uniquely shaped gel pad.

- Silicone pad sits on top of the patellar tendon to relieve pain from patellar tendonitis, patellofemoral syndrome, or Osgood-Schlatter Disease
- Lined with a breathable, hypoallergenic microfleece
- Innovative two-strap design

PRODUCT	PART #
Q Baby Patellar Tendon Strap	5165*

SPINE PRODUCTS


Product Shown: Baja 627 Back Brace





SPINAL BRACING SIZE CHART

- Call: 800-541-2455
- Email: orders@bioskin.com
- Online: www.pro.bioskin.com

SIZING CHART				
	Waist measurement			
	Size	*	Inches	Centimeters
	BAJA™			
	S	1	26 – 30	66 – 76
	M	2	30 – 36	76 – 91
	L	3	36 – 44	91 – 112
	XL	4	44 – 54	112 – 137
	UNIV	4U	26 – 54	66 – 137
	BACK FLEX & BACK PLUS			
	1	1	28 – 34	71 – 86
	3	3	34 – 44	86 – 112
	5	5	44 – 56	112 – 142
	Extension Panels available. Part #97363			
	LUMBAR SUPPORT BRACE			
	XS	0	24 – 26	61 – 66
	S	1	26 – 30	66 – 76
	M	2	30 – 34	76 – 86
	L	3	34 – 38	86 – 97
	XL	4	38 – 42	97 – 107
	XXL	5	42 – 46	107 – 117
Extended Sizes	3XL	6	46 - 50	117 - 127
	4XL	7	50 - 54	127 - 137
	5XL	8	54 - 58	137 - 147

Please contact us for any custom sizing needs.

Product#	Size#	To indicate size, replace * with corresponding number from the sizing chart.
XXXX	*	

BAJA SPINAL BRACES

Baja spinal braces use patented technology to deliver targeted compression for immediate relief from a wide variety of conditions causing low back pain.

- The bilateral pulley system draws the lumbar panel against the soft tissues that support the spine to stabilize the vertebrae in arthritic joints and limit painful micromovements.
- Belt provides circumferential compression to increase intracavitary pressure and relieve pressure on discs.
- Pulleys make it easy to apply firm compression even for people with weak hands or limited dexterity.
- Lightweight, breathable fabrics keep users cool and comfortable.
- Rigid anterior panel
- Each brace includes two hot/cold gel pads that slide into the lumbar pad.
- Available in S-XL or in a universal, trim-to-fit size.



lightweight



dynamic stretch



breathable




hypoallergenic



Product Shown: Baja Back Brace

BAJA SPINAL BRACES DELIVER RESULTS

Personalize the Baja to meet the needs of your patients and your business.



HEAT MOLDABLE LUMBAR PANEL

Made of a polycarbonate-infused ABS nylon, the lumbar panel is heat moldable to an individual patient's anatomy for a customized fit with additional support. When properly customized, the Baja may be billed with Medicare's prefabricated, custom-fitted codes.



CORD COVERS

For a streamlined look, choose the cord cover option when you order Baja products. The cord covers fully enclose the cinching mechanism and improve patient comfort.



CUSTOM LOGOS

Contact our Customer Care team to learn about our custom logo program to promote your business.


CLINICALLY VALIDATED PAIN RELIEF

A study using the bilateral cinching mechanism of the Baja showed an average reduction in VAS scores of 40% and a 41% improvement in Roland Morris Disability Questionnaire scores in patients with chronic low back pain.*

* "Prospective Study of Lumbar Back Brace In an Interventional Pain Practice." Rosenthal, Richard M.D. & Spencer, Shawn, BS. Practical Pain Management, (July/August 2011)


WHAT IS TARGETED COMPRESSION?

Posterior Compression



BioSkin's patented design delivers targeted posterior compression directly to the spine and the soft tissues that support it. When the cinching mechanism is engaged, the lumbar panel is drawn forward, focusing the compression where it is needed for pain relief and support.

Lateral Compression



Other cinch-style back braces have a posterior lacing system with the laces positioned across the spine. When a patient engages this style of brace, the two sides pull together and the compression is focused on the sides of the torso without providing much additional support or compression along the spine.

LEVEL 1

BAJA 627 L0642 OR L0627

Targeted lumbar compression for low-back pain relief.

- Lightweight lumbar support



PRODUCT	PART #
Baja 627	3750*_CC
Baja 627 (without cord covers)	3750*
Baja 627 Universal	37504UCC
Baja 627 RAP	67504URAP

LEVEL 2

BAJA 631 L0648 OR L0631

Targeted compression and pain relief with additional stabilization for the thoracic region of the spine.

- Lumbar and thoracic support



PRODUCT	PART #
Baja 631	3770*_CC
Baja 631 (without cord covers)	3770*
Baja 631 Universal	37704UCC
Baja 631 RAP	67704URAP

LEVEL 3

BAJA 637 L0650 OR L0637

Targeted compression and pain relief with lumbar-thoracic support and sagittal-coronal control.

- Rigid lumbar and thoracic support
- Removable rigid side panels



PRODUCT	PART #
Baja 637	3760*_CC
Baja 637 (without cord covers)	3760*
Baja 637 Universal	37604UCC
Baja 637 RAP	37764URAP



BACK FLEX L0642

BioSkin's original lumbar support with a bilateral cinching mechanism for targeted compression.

- Flexible back panel conforms to the natural curve of the low back
- Rigid front panel for additional stability
- Easily integrate hot/cold therapy into the brace
- Includes two hot/cold gel pads

PRODUCT	PART #
Back Flex	3730*



BACK PLUS L0648

A flexible and supportive back brace with shoulder straps for additional stabilization.

- Bilateral pulley system creates vectored compression to the spine
- Flexible back panel conforms to the natural curve of the low back
- Rigid front panel for postural stability and abdominal control
- Easily integrate hot/cold therapy into the brace
- Taller back panel than the Back Flex
- Shoulder straps included for added stabilization
- Includes two hot/cold gel pads

PRODUCT	PART #
Back Plus	3740*



LUMBAR SUPPORT BRACE A4467

An elastic lumbar wrap.

- Available with a foam oval lumbar pad
- Narrow cut in the front allows and encourages mobility
- Flexible side stays for light stabilization without impeding range of motion
- Vented side panels for thermal regulation

PRODUCT	PART #
Lumbar Support Brace w/ Oval Lumbar Pad	4720*



FOOT & ANKLE PRODUCTS

Product Shown: TriLok®





TriLok Ankle Brace

FOOT AND ANKLE SIZE CHART

Contact Us to Place an Order

- Call: 800-541-2455
- Email: orders@bioskin.com
- Online: www.pro.bioskin.com

	By shoe size					
	TRILOK®					
	Size	*	U.S. Men's	U.S. Women's	Europe	UK
	XS	0		4 – 5.5	33 – 37	2 – 3.5
	S	1	5 – 6.5	6 – 7.5	38 – 39	4 – 5.5
	M	2	7 – 9.5	8 – 11	40 – 42	6 – 8.5
	L	3	10 – 12.5	11.5 – 14	43 – 45	9 – 11.5
	XL	4	13 – 16		46 – 49	12 – 16
	HALLUX CONTROL STRAP					
	Size	*	U.S. Men's	U.S. Women's	Europe	UK
	S	1	5 – 10	6 – 11	38 – 43	4 – 9
	L	3	10.5 – 16	11.5+	44 – 50	9.5 – 16
	WALKING BOOTS					
	Size	*	U.S. Men's	U.S. Women's		
	XS	0	2 – 4	3 – 6		
	S	1	4 – 7	6 – 8		
	M	2	7.5 – 10.5	8 – 11		
	L	3	10.5 – 12.5	11 – 13		
	XL	4	12.5+	13+		
	NIGHT SPLINT					
	Size	*	U.S. Men's	U.S. Women's		
	S - M	2	5 – 9	6 – 10		
	L - XL	4	9 – 14	10 – 15		
	Circumference from heel to top of the foot – back to heel					
	AFTR PRODUCTS					
	Size	*	Inches		Centimeters	
	XS - S	1	11 - 13		27.9 - 33	
	M - L	3	12.5 - 14		31.8 - 35.6	
	XL - XXL	5	13.5 - 15		34.3 - 38.1	

Please contact us for any custom sizing needs.

Product#	Size#	To indicate size, replace * with corresponding number from the sizing chart.
XXXX	*	

TRILOK[®]

L1902

A versatile ankle brace that relieves pain from PTTD and plantar fasciitis, and protects against ankle sprains and peroneal tendon damage.

- FootLok[®] strap supports arch from above
- Inner sleeve provides compression for swelling control and pain relief
- Figure-8 stirrup strap provides hindfoot stabilization
- No laces for easy application
- Soft, antimicrobial fleece against the skin for maximum comfort

PRODUCT	PART #
TriLok	5360*

1. SUBTALAR JOINT CONTROL

- Valgus & Varus Correction
- Ankle Stabilization

2. LATERAL COLUMN UNLOADING

- Ankle Stabilization
- Peroneal Tendon Dysfunction

3. MEDIAL COLUMN UNLOADING

- Posterior Tibial Tendon Dysfunction
- Plantar Fasciitis



DESIGN INSPIRED BY ANATOMY

The patented FootLok Strap wraps around the midfoot and forefoot and anchors to the medial or lateral side of the foot to control motion. This non-stretch strap acts as an external ligament to provide suspension and support to weak and injured tendons and ligaments.



LATERAL APPLICATION

Ankle Sprains &
Peroneal Tendinitis/Tendinosis

MEDIAL APPLICATION

PTTD & Plantar Fasciitis



STO™ L-1902

A simplified Teno-Suspension system provides support for ligament, tendon, and plantar fascia injuries.

- Designed for rapid donning and all-day comfort
- Lateral/Medial ligament protection
- Comfortable, lace free, ankle sprain and roll protection

PRODUCT	PART #
STO	5320*



AFTR DC (DORSAL CLOSURE) L1902

A highly adjustable ankle brace for severe swelling and stabilization after surgery or trauma.

- Wraparound design for easy application
- Adjustable as foot volume decreases
- J-shaped gel pads control swelling around the malleoli
- Gel pads may be chilled for cold therapy
- Trimmable for custom fit
- Non-elastic stirrup straps control hindfoot and limit inversion and eversion (not shown)
- Also available with elastic compression wrap

PRODUCT	PART #
AFTR DC	5386*
AFTR DC with Compression Wrap (Special order, lead times vary)	5386*CW



AFTR (PULL-ON) L1902

A compressive ankle brace for mild to moderate swelling and stabilization after injury.

- Reduces swelling and stabilizes joint after injury
- Trimmable for custom fit
- Non-elastic stirrup straps control hindfoot and limit inversion and eversion

PRODUCT	PART #
AFTR without Gel	5380*
AFTR with Gel	5383*



ANKLE SLEEVE WITH COMPRESSION WRAP

A compressive ankle sleeve with a figure-eight compression wrap for added stability.

- Medical-grade compression
- Lined with breathable microfleece
- Optional compression wrap

PRODUCT	PART #
Ankle Sleeve with Compression Wrap (Special order, lead times vary)	5351*



CMC WRAP 10PK CPT 29540

Provides an attachment point for the Weil Digit Strapping System without covering the ankle.

- Use with the Weil Osteotomy Strap or Hallux Control Strap
- Universal sizing

PRODUCT	PART #
CMC Wrap 10-Pack	95362



WEIL OSTEOTOMY STRAP 10PK CPT 29550

Prevents dorsal drift (floating toe) after Weil Osteotomy by holding the toe in a neutral or plantarflexed position.

- Promotes proper toe alignment
- Use with AFTR braces or the CMC Wrap (pictured)
- Universal sizing

PRODUCT	PART #
Weil Osteotomy Strap 10-Pack	95372



HALLUX CONTROL STRAP 10 PK CPT 29550

Maintains the hallux in a neutral position during critical, soft-tissue healing after bunionectomy.

- Stretches extensor tendons after primary healing to improve range of motion and allow normal gait
- Provides resistance for patients to strengthen flexor hallucis longus for strong toe-off
- Use with AFTR braces or the CMC Wrap (pictured)

PRODUCT	PART #
Hallux Control Strap 10-Pack (Small)	95341T
Hallux Control Strap 10-Pack (Large)	95343T



Product Shown: Aftr

www.bioskin.com



UPPER BODY PRODUCTS


Product Shown: DP2





SIZE CHARTS

Contact Us to Place an Order

- Call: 800-541-2455
- Email: orders@bioskin.com
- Online: www.pro.bioskin.com

	WRIST			
	Circumference of extended wrist			
	Size	*	Inches	Centimeters
	XS – S	1	4 – 6	10 – 15
	M – L	3	6 – 8	15 – 20
	XL – 2XL	5	8 – 10	21 – 25

	THUMB SPICA			
	Length of thumb			
	Size	*	Inches	Centimeters
	XS – M	2	3.25 – 4.25	8 – 11
	L – 2XL	5	4.25 – 5.25	11 – 14

	ELBOW			
	Circumference of largest part of forearm			
	Size	*	Inches	Centimeters
	XS	0	8 – 9	20 – 23
	S	1	9 – 10	23 – 26
	M	2	10 – 11	26 – 28
	L	3	11 – 12	28 – 31
	XL	4	12 – 13	31 – 33
	2XL	5	13 – 14	33 – 36

Please contact us for any custom sizing needs.

Product#	Size#	To indicate size, replace * with corresponding number from the sizing chart.
XXXX	*	



DP2 or DP3 L3908

A cock-up wrist brace made with comfortable materials and adjustable stays.

- Customizable aluminum stays can be bent for comfort and specific positioning or removed if not needed
- Soft microfiber fabric wicks moisture
- Can be applied using only one hand
- Allows full use of the fingers while limiting motion in the wrist

PRODUCT

PART

DP2	5423*L or 5423*R
DP3	5424*L or 5424*R



WRIST THUMB SPICA L3807

Our most supportive wrist brace controls motion in the wrist as well as the CMC and MP joints.

- Allows full use of the fingers while limiting motion in the wrist and thumb joint
- Customizable aluminum stays in the wrist and thumb can be bent for comfort or specific positioning. The palmar stay may be removed if not needed
- Soft microfiber fabric wicks moisture
- Can be applied using only one hand

PRODUCT

PART

Wrist Thumb Spica	5452*L or 5452*R
-------------------	------------------



THUMB SPICA L3923

Low-profile brace controls motion in the CMC and MP joints.

- Customizable aluminum stay bends to fit more comfortably or reposition the thumb
- Soft microfiber fabric wicks moisture during use
- Can be applied using only one hand
- Adjustable compression wrap around the wrist

PRODUCT

PART

Thumb Spica	5450*
-------------	-------



BOOMERANG

Comfortable compression wrist wrap.

- Universal sizing
- Fits either left or right wrist
- Comfortable compression for minor wrist injuries
- Soft microfiber fabric wicks moisture
- Can be applied using only one hand

PRODUCT

PART

Boomerang	54200
-----------	-------



CONCORDE SHOULDER IMMOBILIZER L3960

Adjustable shoulder brace for recovery after injury or surgery requiring diagnosis-specific positioning.

- Abduction bar can be set from neutral to 90 degrees
- External rotation bar can be set from neutral to 90 degrees
- Universal size, moldable waistband and adjustable sling length
- Hypoallergenic, breathable fabric
- Easy application buckle
- Fits left or right arm
- Adjustable height

PRODUCT

PART

Concorde Shoulder Immobilizer

38100



SHOULDER SLING WITH ABDUCTION PILLOW L3960

An adjustable shoulder sling that can be used with or without abduction control.

- Universal size with adjustable sling length
- Use with or without the abduction pillow
- Fits left or right arm
- Includes an optional tension-relieving strap
- Hypoallergenic, breathable fabric
- Easy application buckle for patient compliance

PRODUCT

PART

Shoulder Sling with Abduction Pillow

38200



ELBOW SLEEVE

A compression sleeve for the elbow and forearm.

- Constructed with BioSkin's soft, moisture-wicking materials
- Controls joint swelling and facilitates recovery from minor injuries

PRODUCT

PART

Elbow Sleeve

4210*



TENNIS ELBOW BAND

A low-profile compression band with a gel pad for relief from Tennis Elbow or Golfer's Elbow.

- Constructed with BioSkin's soft, moisture-wicking materials
- Controls joint swelling and facilitates recovery from minor injuries
- Silicone gel pad with strap provides targeted compression to the injured area
- Fits below the crease of the elbow

PRODUCT

PART

Tennis Elbow Band

4250*

NEW! BIOSKIN WRAP

Experience the future of injury management with BioSkin's innovative wrap solution.



WRAP INJURED JOINTS



SECURE A SPLINT



BANDAGE AND IMMOBILIZE



BIOSKIN WRAP

A unique self-adhering, reusable athletic wrap that won't migrate or slip

- Adaptable to any need from small joints to larger areas
- Can be used in conjunction with apparel to wrap areas prone to migration
- Soft-touch micro hook and loop fabric system means no additional fasteners are required
- Creates cohesive bracing solution once applied
- Patent pending calibrated pressure application systems

PRODUCT

PART

Micro 5 - 4" Wrap	PN39054x6
Micro 10 - 3" Wrap	PN39103x6
Micro 10 - 4" Wrap	PN39104x6
Micro 20 - 2" Wrap	PN39202x6
Micro 20 - 4" Wrap	PN39204x6



HIP, THIGH & CALF PRODUCTS




Product Shown: Compression Shorts





SIZE CHARTS

Contact Us to Place an Order

- Call: 800-541-2455
- Email: orders@bioskin.com
- Online: www.pro.bioskin.com

COMPRESSION SHORTS				
	Circumference between hip and top of pelvis			
	Size	*	Inches	Centimeters
	S	1	28 – 32	71 – 81
	M	2	32 – 36	81 – 92
	L	3	36 – 40	92 – 102
	XL	4	40 – 44	102 – 112
	2XL	5	44 – 48	112 – 122

THIGH				
	Mid-thigh circumference			
	Size	*	Inches	Centimeters
	XS	0	16 – 18	41 – 46
	S	1	18 – 20	46 – 51
	M	2	20 – 22	51 – 56
	L	3	22 – 25	56 – 64
	XL	4	25 – 28	64 – 71
	2XL	5	28 – 32	71 – 81
	3XL	6	32 – 36	81 – 91

CALF SLEEVE				
	Largest part of calf circumference			
	Size	*	Inches	Centimeters
	XS	0	12 – 13.5	30 – 34
	S	1	13.5 – 15	34 – 38
	M	2	15 – 16.5	38 – 42
	L	3	16.5 – 18.5	42 – 47
	XL	4	18.5 – 21	47 – 53
	2XL	5	21 – 24	53 – 61
	3XL	6	24 – 27	61 – 69
	4XL	7	27 – 30	69 – 76
	5XL	8	30 – 33	76 – 84

Please contact us for any custom sizing needs.

Product#	Size#	To indicate size, replace * with corresponding number from the sizing chart.
XXXX	*	



ULTIMA COMPRESSION SHORTS

Designed to enhance recovery and improve performance. Optional groin wrap available.

- Hypoallergenic materials
- Optional groin wrap for targeted compression
- Medical-grade compression
- Elastic waistband

PRODUCT	PART #
Compression Shorts	4620*
Compression Shorts with Groin Wrap	4621*



ULTIMA THIGH SLEEVE

Thigh sleeve uses medical-grade compression to relieve pain and enhance recovery.

- Hypoallergenic materials
- SkinLok feature locks sleeve in place
- Medical-grade compression

PRODUCT	PART #
Ultima Thigh Sleeve	4610*



ULTIMA CALF SLEEVE

Designed to enhance recovery and improve performance with medical-grade compression.

- Hypoallergenic materials
- SkinLok feature locks sleeve in place
- Medical-grade compression

PRODUCT	PART #
Ultima Calf Sleeve - Single	4630*



CUSTOM LOGO PROGRAM

Program Details and Requirements

Elevate your brand with the BioSkin Custom Logo Program, a unique opportunity for our valued customers to showcase their identity alongside ours on our innovative medical bracing solutions. This program empowers healthcare professionals, clinics, and organizations to strengthen their presence in patient care by featuring both the BioSkin logo and your custom logo on our high-quality braces, fostering trust and recognition while supporting recovery and mobility. Join us in this collaborative journey to enhance your practice and deliver personalized care with BioSkin's trusted products.

Co-Branding Design

Both the BioSkin logo and your custom logo will be prominently displayed on the braces, ensuring your brand stands out alongside ours without integration into a single design.

Initial Order Lead Time

Allow 4-5 weeks for the first order to be processed and delivered to your location, ensuring meticulous attention to detail.

Artwork Submission

Logos must be submitted as vector files (e.g., .AI, .EPS, or .SVG) to guarantee high-quality reproduction on our products.

Color Guidelines

Custom logos may include color, but black is not permitted as it may not be visible against the fabric; please choose contrasting shades for optimal results.

Size and Eligibility

Heat transfers for custom logos are limited to a maximum of 1.5 inches. Once established, this feature is available for orders totaling \$250 or more.

Cost Structure

A fee of \$3.00 per brace with a custom logo will be applied, reflecting the added value of personalized branding.

Production and Shipping

After the initial setup, logo orders will be assembled within 5-7 business days, with an additional 2-5 business days for standard shipping.

Eligible Products

The following BioSkin bracing products are available for custom logo application:

- Trilok Ankle Brace
- AFTR Ankle Braces
- Tennis Elbow Band
- Concord Shoulder
- Elbow Sleeve
- Thumb Spica
- Wrist/Thumb Spica
- DP2/DP3
- All sewn Knee Braces
- Baja Back Braces with Cord Cover
- Shoulder sling with Abduction Pillow





Pro.BioSkin.com

CALL: 800-541-2455 | FAX: 541-482-2119 | EMAIL: orders@bioskin.com



BIOSKIN™