



BioSkin™



# INNOVATIVE BRACING SOLUTIONS

Please click any of the links below to jump to that section.

- [Knee](#)
- [Spine](#)
- [Foot & Ankle](#)
- [Upper Extremity](#)
- [Lower Extremity](#)





*TriLok™ Ankle Brace*

[pro.bioskin.com](http://pro.bioskin.com)

CALL: 800-541-2455 | FAX: 541-482-2119 | EMAIL: [orders@bioskin.com](mailto:orders@bioskin.com)

# 25 YEARS OF INNOVATION

BioSkin is both a family-owned company and a world leader in the design of orthopedic braces and sports medicine products. We've grown our business on the development of advanced materials and sophisticated bracing solutions that decrease pain, restore mobility, and improve lives. We remain committed to helping people overcome challenges and stay active in every stage of life.

## INTERNATIONAL

We have distributors located all over the world.



# THE BIOSKIN DIFFERENCE

Twenty-five years ago, we created a breathable, lightweight, and compressive alternative to the neoprene traditionally used in braces and called it BioSkin.

Today, we select from a family of proprietary materials with varied technical characteristics when developing new products. These materials allow us to fine-tune every aspect of the design to maximize the effectiveness and comfort of each brace.



lightweight



dynamic stretch



breathable



hypoallergenic

## BIOSKIN

- Cool
- Breathable
- Hypoallergenic
- 1 millimeter

## NEOPRENE


- Hot
- Nonbreathable
- Common skin irritant
- 7+ millimeters





# DESIGN INSPIRED BY ANATOMY

We design our braces to work in harmony with the human body, taking our cues from anatomical structures and biomechanical function. The result is products that encourage healthy movement and increase the quality of life for people who rely on them.



# A PARTNERSHIP FOR EXCEPTIONAL CARE

We know your goal is to provide exceptional care to your patients. Here are some of the ways our products help you do that:

- Thoughtfully-designed products prevent injury, reduce pain, and restore mobility
- Hypoallergenic materials are dry, cool, and comfortable for all-day use
- A wide range of sizes and design options provide the best fit possible
- Quality components are lightweight and durable for comfortable products that last
- Attention to details like carefully placed seams and non-constricting straps increase patient compliance







[pro.bioskin.com](http://pro.bioskin.com)

CALL: 800-541-2455 | FAX: 541-482-2119 | EMAIL: [orders@bioskin.com](mailto:orders@bioskin.com)



# KNEE BRACES


**PERFORMANCE**  
**PREVENTION**  
**RECOVERY**

The BioSkin logo, featuring a stylized blue and white wave icon above the brand name "BioSkin" in a sans-serif font.The BioSkin logo, featuring a stylized blue and white wave icon above the brand name "BioSkin" in a sans-serif font.

# EASILY ORDER BIOSKIN PRODUCTS

Contact us to begin the ordering process:

- Call: 800-541-2455
- Email: [orders@bioskin.com](mailto:orders@bioskin.com)
- Online: [www.pro.bioskin.com](http://www.pro.bioskin.com)

KNEE PRODUCTS						
	Circumference with knee fully extended					
	Size	*	Inches		Centimeters	
			mid – patella	7" above	mid – patella	18cm above
	XS	0	12 – 13	16 – 18	31 – 33	41 – 46
	S	1	13 – 14	18 – 20	33 – 36	46 – 51
	M	2	14 – 16	20 – 22	36 – 41	51 – 56
	L	3	16 – 17	22 – 25	41 – 43	56 – 64
	XL	4	17 – 19	25 – 28	43 – 46	64 – 71
	XXL	5	19 – 21	28 – 32	46 – 54	71 – 81
	Extended Sizes	3XL	6	21 – 23	32 – 36	54 – 58
4XL		7	23 – 25	36 – 40	58 – 64	91 – 102
5XL		8	25 – 27	40 – 44	64 – 69	102 – 112

Product#	Size#	To indicate size, replace * with corresponding number from the sizing chart.
XXXX	*	

Please contact us for any custom sizing needs.

## PREMIUM MATERIAL OPTIONS

Some of our braces are available in more than one material, making it even easier to fit every patient. Braces made with Ultima offer a higher level of compression, while braces made with Stratus have less compression and a more comfortable fit that is easier to pull on. When available, these options and associated part numbers are noted.

# EVERY DETAIL MATTERS

A brace can't help if the patient won't wear it, so we obsess over details to make sure our products are both supportive and comfortable.

## PROPRIETARY MATERIALS

- Medical-grade compression reduces pain and swelling
- Moisture activates the SkinLok™ feature and prevents migration

## COMFORTABLE STRAPS

- Elastic straps allow muscle expansion and minimize migration
- Die-cut velcro tabs eliminate sharp corners

## BUILT FOR COMFORT

- No seams behind the knee

## HINGES/SUPPORTS

- High quality hinges are both durable and lightweight



lightweight



dynamic stretch



breathable



hypoallergenic



# PERFECT FIT

Every leg is different and one-size-fits-all braces sacrifice efficacy and compliance. That's why we have multiple configurations as well as extended sizing from XS - 5XL.



## WRAPAROUND DESIGN

Many of our braces are available with our unique wraparound design. The entire front panel is Velcro compatible, making it easy to achieve a perfect fit on virtually any leg. This is also a great option for patients with limited flexibility or hand strength.

## PULL-ON DESIGN

Our pull-on knee braces are made with a highly compressive lycra-based fabric to maximize pain relief, control swelling, and enhance proprioception.



# CROSSFIRE™

Any brace with CrossFire in the name incorporates innovative design features to maximize the therapeutic benefits of compression and increase comfort.



## CROSSFIRE HYBRID WRAPAROUND DESIGN

An inner sleeve facilitates donning by holding the brace in place so patients have both hands available to adjust the fit. The sleeve provides compression to the entire knee to control swelling and deliver pain relief while the wraparound component of the brace also provides more sizing flexibility and a perfect fit.



## ULTRA LIGHTWEIGHT MATERIAL BEHIND THE KNEE

CrossFire wraparound and pull-on braces feature our exceptionally lightweight Cirrus material behind the knee. This back panel is sewn in a slightly flexed position, virtually eliminating bunching behind the knee.



## CROSSFIRE PULL-ON DESIGN

Some CrossFire braces are also available in a pull-on version that features the same lightweight and compressive material over the patella as shown in the wraparound version's donning sleeve. These braces also feature a more generous cut above the knee to accommodate a wider variety of leg shapes.

A person is shown from the waist down, wearing a black knee brace on their right knee. The brace has a white logo that says "BioSkin" and a dynamic spiral strap. The person is wearing blue athletic shoes and is standing on a stone step. The background is a blurred green foliage. A large blue curved shape is on the right side of the image.

# ELIMINATE OA KNEE PAIN

Dynamic spiral strap  
gently pulls the joint  
open during extension,  
applying the unloading  
force only when needed.



# CROSSFIRE OA

L1851/L1843

## CROSSFIRE OA - DYNAMIC UNLOADER

**Offloads the medial or lateral compartment to reduce osteoarthritis or meniscus pain. Also addresses patellofemoral pain.**

- Spiral strap works in tandem with the hinge, calf, and thigh to create a three-point leverage system that opens the joint during extension. The strap relaxes during flexion.
- Magnetic buckle attachment on the spiral strap allows for quick and easy application and removal of the brace.
- Optional calf and thigh straps provide additional support.

PRODUCT	PART #
CrossFire OA	5190* (L or R)

*Note: Use "L" for the left side brace and "R" for the right side.*



*Optional patella traction strap can adjust the tension on the patella to correct maltracking and relieve pain.*

## ROM HINGES: KNEE STABILIZATION WITH MOTION CONTROL

These braces provide high-level stabilization and limit flexion and extension with adjustable range-of-motion hinges.



### GLADIATOR™

**Our most supportive hinged knee brace.**

- Non-elastic straps for maximum stabilization
- Rigid, aluminum, 15" uprights may be bent for a customized fit (L1832)
- Gladiator hinge: ROM hinge with flexion/extension stops at 15°, 30°, 60°, and 90°
- Pull-on foam buttress around the patella for additional stabilization

PRODUCT	PART #
Gladiator   Wraparound	5136*
Gladiator   Pull-on (Ultima)	4134*
Gladiator   Pull-on (Stratus)	6136*



### CROSSFIRE ROM

**Exceptionally comfortable with a lighter weight hinge.**

- Rigid, aluminum uprights are flared to contour to the leg's natural shape and may be bent for additional customization (L1832)
- SlimLok hinge: Lightweight, adjustable ROM hinge with three flexion and extension stops
- Pull-on version available with a silicone gel pad around the patella to reduce swelling and decrease knee pain

PRODUCT	PART #
CrossFire ROM   Wraparound	5180*
CrossFire ROM   Pull-on	6180*
CrossFire ROM   Pull-on with Gel	6187*



### Q LOK™ ROM

**Patellar control in a highly supportive hinged knee brace.**

- Traction strap controls patellar glide and tilt. Fits left or right leg
- Non-elastic straps control medial/lateral and anterior/posterior movement of the joint
- Rigid, aluminum uprights are flared to contour to the leg's natural shape and can be bent for additional customization (L1832)
- SlimLok hinge: Lightweight, adjustable ROM hinge with three flexion and extension stops at 0°, 15°, and 30°

PRODUCT	PART #
Q Lok ROM   Wraparound (interior donning sleeve)	5166*
Q Lok ROM   Pull-on	6166*

## LIGHTWEIGHT ALUMINUM HINGES: MODERATE STABILIZATION

Bicentric, aluminum hinges are lightweight yet supportive. Made of aircraft-grade aluminum for its unique strength-to-weight ratio.



### CROSSFIRE ALUMINUM

**Exceptionally comfortable premium hinged knee brace.**

**We've updated our top-selling knee brace to include the CrossFire design (p. 5).**

- Inner donning sleeve for ease of application and light compression over the patella
- Exceptionally lightweight material is sewn in tension behind the knee to prevent uncomfortable bunching

PRODUCT	PART #
CrossFire Aluminum   Wraparound	5186*
CrossFire Aluminum   Pull-on	6186*



### HINGED KNEE SKIN™

**Lightweight, hinged knee brace for general joint stabilization.**

- Bicentric, aluminum hinges support ligaments and allow natural motion of the knee
- Pull-on version includes foam buttress around the patella for additional stabilization

PRODUCT	PART #
Hinged Knee Skin   Wraparound	5133*
Hinged Knee Skin   Pull-on (Ultima)	4130*
Hinged Knee Skin   Pull-on (Stratus)	6130*

*Unless noted, all of our braces are universal left or right.*



## CONFORMA™ HINGE: FLEXIBLE, LIGHTWEIGHT SUPPORT

The Conforma hinge is made of lightweight, flexible nylon that conforms to the shape of any leg. Surprisingly lightweight, comfortable, and durable.



### Q BRACE™ WITH CONFORMA HINGE

**A patellar tracking brace with multidirectional patella strap and lightweight hinges.**

- Universal design (Q strap works on either knee and can track in multiple directions)
- Prevents subluxation
- May also be worn as a hinged knee brace without the Q strap
- Also available without a hinge (A4467)

PRODUCT	PART #
Q Brace w/ Conforma Hinge   Wraparound	5138*
Q Brace w/ Conforma Hinge   Pull-on	4136*



### PREMIUM J™ WITH CONFORMA HINGE

**A simple and comfortable solution for patellofemoral pain and joint stabilization.**

- J-shaped foam buttress creates a lateral barrier for the patella to promote proper alignment, prevent subluxation, and provide pain relief
- Elastic patellar control straps and closed patella design provide lateral pressure and anterior compression to increase surface contact in the patellofemoral joint
- Also available without a hinge (A4467)

PRODUCT	PART #
Premium J w/ Conforma Hinge	6138*(L or R)

*Note: over. Use "L" for the left side brace and "R" for the right side.*



### PATELLA STABILIZER WITH CONFORMA HINGE

**A lightweight hinged knee brace with a foam buttress around the patella for joint stabilization.**

- Foam ring around the patella for stabilization
- Open patella design
- Also available without a hinge (A4467)

PRODUCT	PART #
Patella Stabilizer w/ Conforma Hinge (Ultima)	4132*
Patella Stabilizer w/ Conforma Hinge (Stratus)	6132*



## CROSSFIRE CONFORMA

**A lightweight hinged knee brace for pain relief and light stabilization.**

- Lightweight, flexible hinges provide sagittal control without sacrificing comfort
- Now also available with silicone gel ring around the patella to control swelling and maximize pain relief



PRODUCT	PART #
CrossFire Conforma   Wraparound	5185*
CrossFire Conforma   Pull-on	6185*
CrossFire Conforma   Pull-on with Gel	6188*

## VISCO KNEE SKIN™ WITH CONFORMA HINGE

**An all-purpose knee brace with a silicone gel ring perfect for relieving general knee pain and swelling.**

- Soft, silicone gel ring around the patella reduces swelling and increases compression around the joint
- Also available without a hinge (A4467)



PRODUCT	PART #
Visco Knee Skin w/ Conforma Hinge	6137*

## NON-HINGED KNEE SUPPORTS

Medical-grade compression for pain relief



### TEN-7™

**A targeted compression sleeve, uniquely designed to relieve chronic knee pain.**

- Provides targeted compression to the anatomical structures proximal to the knee
- Worn comfortably and discreetly on the thigh
- Relieves pain from osteoarthritis, meniscus damage and other chronic conditions



### KNEE SKIN™

**A simple compression sleeve for comfortable support and pain relief.**

- Flexible stays on each side provide light stabilization to the knee
- Compression around the knee reduces pain and swelling
- Available with or without stabilization straps
- Available in an open or closed patella design

PRODUCT	PART #
TEN-7	5146*

PRODUCT	PART #
Knee Skin (without straps)	4110* (O OR C)
Knee Skin (with straps)	4111* (O OR C)

*Note: Use "O" for open patella and "C" for closed patella option*



### UNDER SKIN™

**A breathable, compression sleeve for use under a rigid brace to prevent migration and protect skin.**

- Constructed with lightweight yet compressive fabric for maximum comfort and enhanced recovery



### Q BABY™

**A premium patellar tendon strap with a uniquely shaped gel pad.**


- Silicone pad sits on top of the patellar tendon to relieve pain from patellar tendonitis, patellofemoral syndrome, or Osgood-Schlatter Disease
- Lined with a breathable, hypoallergenic microfleece
- Innovative two-strap design

PRODUCT	PART #
Under Skin	4140*

PRODUCT	PART #
Q Baby Patellar Tendon Strap	5165*




# SPINE PRODUCTS



# THREE WAYS TO ORDER

- Call: 800-541-2455
- Email: [orders@bioskin.com](mailto:orders@bioskin.com)
- Online: [www.pro.bioskin.com](http://www.pro.bioskin.com)

SIZING CHART				
	Waist measurement			
	Size	*	Inches	Centimeters
	BAJA™			
	S	1	26 – 30	66 – 76
	M	2	30 – 36	76 – 91
	L	3	36 – 44	91 – 112
	XL	4	44 – 54	112 – 137
	UNIV	4U	26 – 54	66 – 137
	BACK FLEX & BACK PLUS			
	1	1	28 – 34	71 – 86
	3	3	34 – 44	86 – 112
	5	5	44 – 56	112 – 142
	Extension Panels available. Part #97363			
	BACK SKIN			
	XS	0	24 – 26	61 – 66
	S	1	26 – 30	66 – 76
	M	2	30 – 34	76 – 86
	L	3	34 – 38	86 – 97
	XL	4	38 – 42	97 – 107
	XXL	5	42 – 46	107 – 117
Extended Sizes	3XL	6	46 - 50	117 - 127
	4XL	7	50 - 54	127 - 137
	5XL	8	54 - 58	137 - 147

Product#	Size#	To indicate size, replace * with corresponding number from the sizing chart.
XXXX	*	

Please contact us for any custom sizing needs.



# BAJA SPINAL BRACES

Baja spinal braces use patented technology to deliver targeted compression for immediate relief from a wide variety of conditions causing low back pain.

- The bilateral pulley system draws the lumbar panel against the soft tissues that support the spine to stabilize the vertebrae in arthritic joints and limit painful micromovements.
- Belt provides circumferential compression to increase intracavitary pressure and relieve pressure on discs.
- Pulleys make it easy to apply firm compression even for people with weak hands or limited dexterity.
- Lightweight, breathable fabrics keep users cool and comfortable.
- Each brace includes two hot/cold gel pads that slide into the lumbar pad.
- Available in S-XL or in a universal, trim-to-fit size.



lightweight



dynamic stretch



breathable



hypoallergenic

*Baja 627  
(See page 5)*



# REDESIGNED FOR IMPROVED PERFORMANCE

The redesigned Baja line incorporates three new features that increase comfort, efficacy and improve the overall look.



## HEAT MOLDABLE LUMBAR PANEL

Now made of a polycarbonate-infused ABS nylon, the new lumbar panel is heat moldable to an individual patient's anatomy for a customized fit with additional support. When properly customized, the Baja may be billed with Medicare's prefabricated, custom-fitted codes.



## CORD COVERS

The redesigned Baja line fully encloses the cords of the cinching mechanism for improved patient comfort and a more streamlined look.



## CUSTOM LOGOS

Order 20 units or more and add your logo to the brace to promote your business.


# CLINICALLY VALIDATED PAIN RELIEF

A study using the bilateral cinching mechanism of the Baja showed an average reduction in VAS scores of 40% and a 41% improvement in Roland Morris Disability Questionnaire scores in patients with chronic low back pain.\*

\* "Prospective Study of Lumbar Back Brace In an Interventional Pain Practice. " Rosenthal, Richard M.D. & Spencer, Shawn, BS. Practical Pain Management, (July/August 2011)


## WHAT IS TARGETED COMPRESSION?

### Posterior Compression



BioSkin's patented design delivers targeted posterior compression directly to the spine and the soft tissues that support it. When the cinching mechanism is engaged, the lumbar panel is drawn forward, focusing the compression where it is needed for pain relief and support.

### Lateral Compression



Other cinch-style back braces have a posterior lacing system with the laces positioned across the spine. When a patient engages this style of brace, the two sides pull together and the compression is focused on the sides of the torso without providing much additional support or compression along the spine.



## LEVEL 1

### BAJA 627

Targeted lumbar compression for low-back pain relief.

- Lightweight lumbar support
- PDAC-approved for **L0642** or **L0627**

#### PRODUCT

#### PART #

Baja 627

3750\*\_CC

Baja 627 (without cord covers)

3750\*



## LEVEL 2

### BAJA 631

Targeted compression and pain relief with additional stabilization for the thoracic region of the spine.

- Lumbar and thoracic support
- PDAC-approved for **L0648** or **L0631**

#### PRODUCT

#### PART #

Baja 631

3770\*\_CC

Baja 631 (without cord covers)

3770\*



## LEVEL 3

### BAJA 637

Targeted compression and pain relief with lumbar-thoracic support and sagittal-coronal control.

- Rigid lumbar and thoracic support
- Removable rigid side panels
- PDAC approved for **L0650** or **L0637**

#### PRODUCT

#### PART #

Baja 637

3760\*\_CC

Baja 637 (without cord covers)

3760\*





## BACK FLEX L0642

**BioSkin's original lumbar support with a bilateral cinching mechanism for targeted compression.**

- Flexible back panel conforms to the natural curve of the low back
- Rigid front panel for additional stability
- Easily integrate hot/cold therapy into the brace
- Two hot/cold gel pads included with purchase

### PRODUCT

### PART #

Back Flex

3730\*



## BACK PLUS L0648

**A flexible and supportive back brace with shoulder straps for additional stabilization.**

- Bilateral pulley system creates vectored compression to the spine
- Flexible back panel conforms to the natural curve of the low back
- Rigid front panel for postural stability and abdominal control
- Easily integrate hot/cold therapy into the brace
- Gel pads can be heated in the microwave or cooled in the freezer
- Taller back panel than the Back Flex
- Shoulder straps included for added stabilization
- Two hot/cold gel pads included with purchase

### PRODUCT

### PART #

Back Plus

3740\*



## BACK SKIN

**An elastic lumbar wrap with additional support options.**

- Available with foam lumbar cushion or flexible lumbar pad
- Narrow cut in the front allows and encourages mobility
- Flexible side stays for light stabilization without impeding range of motion
- Vented side panels for thermal regulation

### PRODUCT

### PART #

Back Skin (HCPCS: A4467)

4710\*

Back Skin w/ Foam Lumbar Pad (HCPCS: A4467)

4720\*

Back Skin w/ Flexible Lumbar Support (HCPCS: A4467)

4734\*

Back Skin w/ Foam Pad & Flexible Support (HCPCS: A4467)

4724\*


# FOOT & ANKLE PRODUCTS




# EASILY ORDER BIOSKIN PRODUCTS

Contact us to begin the ordering process:

- Call: 800-541-2455
- Email: [orders@bioskin.com](mailto:orders@bioskin.com)
- Online: [www.pro.bioskin.com](http://www.pro.bioskin.com)



By shoe size					
TRILOK™					
Size	*	U.S. Men's	U.S. Women's	Europe	UK
XS	0		4 – 5.5	35 – 37	2 – 3.5
S	1	5 – 6.5	6 – 7.5	38 – 39	4 – 5.5
M	2	7 – 9.5	8 – 11	40 – 42	6 – 8.5
L	3	10 – 12.5	11.5 – 14	43 – 45	9 – 11.5
XL	4	13 – 16		46 – 49	12 – 16
HCS - HALLUX CONTROL STRAP					
Size	*	U.S. Men's	U.S. Women's	Europe	UK
S	1	5 – 10	6 – 11	38 – 43	4 – 9
L	3	10.5 – 16	11.5+	44 – 50	9.5 – 16
WALKING BOOTS					
Size	*	U.S. Men's	U.S. Women's		
XS	0	2 – 4	3 – 6		
S	1	4 – 7	6 – 8		
M	2	7.5 – 10.5	8 – 11		
L	3	10.5 – 12.5	11 – 13		
XL	4	12.5+	13+		
NIGHT SPLINT					
Size	*	U.S. Men's	U.S. Women's		
S - M	2	5 – 9	6 – 10		
L - XL	4	9 – 14	10 – 15		



Circumference from heel to top of the foot – back to heel			
AFTR PRODUCTS			
Size	*	Inches	Centimeters
XS - S	1	11 - 13	27.9 - 33
M - L	3	12.5 - 14	31.8 - 35.6
XL - XXL	5	13.5 - 15	34.3 - 38.1



# EVERY DETAIL MATTERS

A brace isn't helpful if the patient won't wear it, so we obsess over details to make sure our products are both supportive and comfortable.



lightweight



dynamic stretch



breathable



hypoallergenic

## INNOVATIVE MATERIALS

- Medical-grade compression reduces pain and swelling

## DESIGN INSPIRED BY ANATOMY

- Mimics the body's natural biomechanics to provide exceptional support
- Limits harmful or painful motion while encouraging mobility

## MADE WITH YOUR PATIENTS IN MIND

- Fits in most shoes
- Easy care and cleaning
- Die-cut velcro tabs eliminate sharp corners
- Adjustable for a customized fit



TriLok Ankle Brace

# TRILOK™

## L1902

A versatile ankle brace that relieves pain from PTTD and plantar fasciitis, and protects against ankle sprains and peroneal tendon damage.

- FootLok™ strap supports arch from above
- Inner sleeve provides compression for swelling control and pain relief
- Figure-8 stirrup strap provides hindfoot stabilization
- No laces for easy application
- Soft, antimicrobial fleece against the skin for maximum comfort

### PRODUCT

TriLok

### PART #

5360\*

### 1. SUBTALAR JOINT CONTROL

- Valgus & Varus Correction
- Ankle Stabilization

### 2. LATERAL COLUMN UNLOADING

- Ankle Stabilization
- Peroneal Tendon Dysfunction

### 3. MEDIAL COLUMN UNLOADING


- Posterior Tibial Tendon Dysfunction
- Plantar Fasciitis





## DESIGN INSPIRED BY ANATOMY

The patented FootLok Strap wraps around the midfoot and forefoot and anchors to the medial or lateral side of the foot to control motion. This non-stretch strap acts as an external ligament to provide suspension and support to weak and injured tendons and ligaments.





## AFTR DC (DORSAL CLOSURE)

**A highly adjustable ankle brace for severe swelling and stabilization after surgery or trauma.**

- Wraparound design for easy application
- Adjustable as foot volume decrease
- J-shaped gel pads control swelling around the malleoli
- Gel pads may be chilled for cold therapy
- Trimmable for custom fit
- Non-elastic stirrup straps control hindfoot and limit inversion and eversion (not shown)
- Also available with elastic compression wrap



PRODUCT	PART #
AFTR DC	5386*
AFTR DC with Compression Wrap	5386*CW

## AFTR (PULL-ON)

**A compressive ankle brace for mild to moderate swelling and stabilization after injury.**

- Reduces swelling and stabilizes joint after injury
- Trimmable for custom fit
- Non-elastic stirrup straps control hindfoot and limit inversion and eversion
- Also available with elastic compression wrap



PRODUCT	PART #
AFTR without Gel	5380*
AFTR with Compression Wrap	5380*CW
AFTR with Gel	5383*
AFTR with Gel and Compression Wrap	5383*CW

## ANKLE SKIN

**A simple compression ankle sleeve for edema control and light stabilization.**

- Pulls on like a sock
- Also available with elastic compression wrap



PRODUCT	PART #
Ankle Skin	5350*
Ankle Skin with Compression Wrap	5351*

**CMC WRAP 10PK**

Provides an attachment point for the Weil Digit Strapping System without covering the ankle.

- Use with the Weil Osteotomy Strap or Hallux Control Strap
- Universal sizing

PRODUCT	PART #
CMC Wrap 10-Pack	95362

**WEIL OSTEOTOMY STRAP 10PK**

Prevents dorsal drift (floating toe) after Weil Osteotomy by holding the toe in a neutral or plantarflexed position.

- Promotes proper toe alignment
- Use with AFTR braces or the CMC Wrap (pictured)
- Universal sizing

PRODUCT	PART #
Weil Osteotomy Strap 10-Pack	95372

**HALLUX CONTROL STRAP 10 PK**

Maintains the hallux in a neutral position during critical, soft-tissue healing after bunionectomy.

- Stretches extensor tendons after primary healing to improve range of motion and allow normal gait
- Provides resistance for patients to strengthen flexor hallucis longus for strong toe-off
- Use with AFTR braces or the CMC Wrap (pictured)

PRODUCT	PART #
HCS 10-Pack (Small)	95341T
HCS 10-Pack (Large)	95343T

## NIGHT SPLINT

A protective splint to promote proper healing of soft tissue during sleep.



- Anti-slip sole
- Adjustable tension strap
- Pneumatic pad for added support
- Articulating hinge design for functional adjustability
- Wide and cushioned metatarsal area for maximal toe comfort

PRODUCT	PART #
Night Splint	3320*

## L4361 & L4387

### PNEUMATIC AND NON-PNEUMATIC WALKING BOOTS - SHORT & TALL

A cushioned walking boot to protect the foot post surgery or traumatic injury.



*Non-Pneumatic Boot*

- Taller boot design increases leverage
- Rounded sole accommodates natural gait
- Shock Pod absorbs over 50% of the force of the heel strike phase
- Anti-slip lip allows for easy elevation
- Soft, anti-microbial fleece against the skin
- Pneumatic cushion for swelling control and pain relief
- Adjustable straps for a customized fit
- Also available in non-pneumatic



*Pneumatic Boot*

PRODUCT	PART #
Non-Pneumatic Boot Short	3311*
Non-Pneumatic Boot Tall	3312*
Pneumatic Boot Short	3303*
Pneumatic Boot Tall	3310*





# UPPER BODY





# HOW TO ORDER

Contact us to begin the ordering process:

- Call: 800-541-2455
- Email: [orders@bioskin.com](mailto:orders@bioskin.com)
- Online: [www.pro.bioskin.com](http://www.pro.bioskin.com)

	ELBOW			
	Circumference of largest part of forearm			
	Size	*	Inches	Centimeters
	XS	0	8 – 9	20 – 23
	S	1	9 – 10	23 – 26
	M	2	10 – 11	26 – 28
	L	3	11 – 12	28 – 31
	XL	4	12 – 13	31 – 33
	2XL	5	13 – 14	33 – 36

	WRIST			
	Circumference of extended wrist			
	Size	*	Inches	Centimeters
	XS – S	1	4 – 6	10 – 15
	M – L	3	6 – 8	15 – 20
	XL – 2XL	5	8 – 10	21 – 25
	THUMB SPICA			
	Length of thumb			
	Size	*	Inches	Centimeters
	XS – M	2	3.25 – 4.25	8 – 11
	L – 2XL	5	4.25 – 5.25	11 – 14

## EVERY DETAIL MATTERS

A brace isn't helpful if the patient won't wear it, so we obsess over details to make sure our products are both supportive and comfortable.

### INNOVATIVE MATERIALS

- Skin-friendly materials are soft, hypoallergenic, and wick moisture



lightweight



dynamic stretch



breathable



hypoallergenic

### ERGONOMIC DESIGNS

- All wrist braces can easily be put on with one hand

### CUSTOMIZED FIT

- Bend the aluminum stays in the wrist braces to the desired position
- Position the silicone gel pad in the elbow sleeves for targeted compression

# L3908

## DP2 or DP3

A cock-up wrist brace made with comfortable materials and adjustable stays.



DP2 - 6 ½"



DP3 - 8 ½"

- Customizable aluminum stays can be bent for comfort and specific positioning or removed if not needed
- Soft microfleece fabric wicks moisture
- Can be applied using only one hand
- Allows full use of the fingers while limiting motion in the wrist

### PRODUCT

### PART #

DP2

5423\*L or 5423\*R

DP3

5424\*L or 5424\*R

# L3807 & L3923

## WRIST THUMB SPICA

Our most supportive wrist brace controls motion in the wrist as well as the CMC and MP joints.



- Customizable aluminum stays in the wrist and thumb can be bent for comfort or specific positioning. The palmar stay may be removed if not needed
- Soft microfleece fabric wicks moisture
- Can be applied using only one hand
- Allows full use of the fingers while limiting motion in the wrist and thumb joint

### PRODUCT

### PART #

Wrist Thumb Spica

5452\*L or 5452\*R

## THUMB SPICA

Low-profile brace controls motion in the CMC and MP joints.



- Customizable aluminum stay bends to fit more comfortably or reposition the thumb
- Soft microfleece fabric wicks moisture during use
- Can be applied using only one hand
- Adjustable compression wrap around the wrist

### PRODUCT

### PART #

Thumb Spica

5450\*L or 5450\*R

# A4467

## BOOMERANG

Comfortable compression wrist wrap.



- Universal sizing
- Fits either left or right wrist
- Comfortable compression for minor wrist injuries
- Soft microfleece fabric wicks moisture
- Can be applied using only one hand

### PRODUCT

### PART #

Boomerang

54200



## ELBOW SKIN

A compression sleeve for the elbow and forearm.

- Constructed with BioSkin's soft, moisture-wicking materials
- Controls joint swelling and facilitates recovery from minor injuries

### PRODUCT

### PART #

Elbow Skin

4210\*

## TENNIS ELBOW SKIN

A compression sleeve with a gel pad for targeted relief from Tennis Elbow or Golfer's Elbow.

- Constructed with BioSkin's soft, moisture-wicking materials
- Controls joint swelling and facilitates recovery from minor injuries
- Silicone gel pad with strap provides targeted compression to the injured area

### PRODUCT

### PART #

Tennis Elbow Skin

4220\*

## TENNIS ELBOW BAND

A low-profile compression band with a gel pad for relief from Tennis Elbow or Golfer's Elbow.

- Constructed with BioSkin's soft, moisture-wicking materials
- Controls joint swelling and facilitates recovery from minor injuries
- Silicone gel pad with strap provides targeted compression to the injured area
- Fits entirely below the crease of the elbow

### PRODUCT

### PART #

Tennis Elbow Band

4250\*

## CONCORDE SHOULDER IMMOBILIZER

Adjustable shoulder brace for recovery after injury or surgery requiring diagnosis-specific positioning.

- Abduction bar can be set from neutral to 90 degrees
- External rotation bar can be set from neutral to 90 degrees
- Universal size, moldable waistband and adjustable sling length
- Hypoallergenic, breathable fabric
- Easy application buckle
- Fits left or right arm
- Adjustable height

### PRODUCT

### PART #

Concorde Shoulder Immobilizer

38100



# HIP, THIGH & CALF PRODUCTS





## COMPRESSION SHORTS

Designed to enhance recovery and improve performance. Optional groin wrap available.

- Hypoallergenic materials
- Optional groin wrap for targeted compression
- Medical-grade compression
- Elastic waistband

PRODUCT	PART #
Compression Shorts	4620*
Compression Shorts with Groin Wrap	4621*



## THIGH SKIN

Thigh sleeve uses medical-grade compression to relieve pain and enhance recovery.

- Hypoallergenic materials
- SkinLok feature locks sleeve in place
- Medical-grade compression
- Available with or without targeted cinch strap

PRODUCT	PART #
Thigh Skin	4610*
Thigh Skin with Cinch Strap	4611*



## CALF SKIN

Designed to enhance recovery and improve performance with medical-grade compression.

- Hypoallergenic materials
- SkinLok feature locks sleeve in place
- Medical-grade compression
- Available as a single sleeve or pair

PRODUCT	PART #
Calf Skin - Single	4630*
Calf Skin - Pair	4630*P

COMPRESSION SHORTS				
Circumference between hip and top of pelvis				
Size	*	Inches	Centimeters	
S	1	28 – 32	71 – 81	
M	2	32 – 36	81 – 92	
L	3	36 – 40	92 – 102	
XL	4	40 – 44	102 – 112	
2XL	5	44 – 48	112 – 122	

THIGH				
Mid-thigh circumference				
Size	*	Inches	Centimeters	
XS	0	16 – 18	41 – 46	
S	1	18 – 20	46 – 51	
M	2	20 – 22	51 – 56	
L	3	22 – 25	56 – 64	
XL	4	25 – 28	64 – 71	
2XL	5	28 – 32	71 – 81	
3XL	6	32 – 36	81 – 91	

CALF				
Largest part of calf circumference				
Size	*	Inches	Centimeters	
S	1	11 – 12	28 – 30	
M	2	12 – 14	30 – 35	
L	3	14 – 16	35 – 41	
XL	4	16 – 18	41 – 46	
2XL	5	18 – 20	46 – 51	
3XL	6	20 – 22	51 – 56	



UPPER BODY



CONCORDE SHOULDER  
IMMOBILIZER

38100 | HCPCS: L3960



BOOMERANG™

54200 | HCPCS: A4467



DP2™ 6 1/2" COCK-UP  
WRIST

5423\*L or R | HCPCS: L3908



DP3™ 8 1/2" COCK-UP  
WRIST

5424\* L or R | HCPCS: L3908



ELBOW SKIN™

4210\* | HCPCS: A4467



TENNIS ELBOW BAND™

4250\* | HCPCS: A4467



TENNIS ELBOW SKIN™

4220\* | HCPCS: A4467



THUMB SPICA


5450\* | HCPCS: L3923



WRIST/THUMB SPICA

5452\* L or R | HCPCS: L3807

PRODUCT  
REFERENCE  
GUIDE



pro.bioskin.com  
CALL: 800-541-2455  
EMAIL: assist@bioskin.com  
©2020 BioSkin® 11/2020 PN: 12659B



FOOT & ANKLE

AFTR™



5380\* | HCPCS: L1902

AFTR DC™



5386\* | HCPCS: L1902

CALIBRATED MIDFOOT  
COMPRESSION WRAP™



\*shown with Weil  
Osteotomy Strap  
(sold separately)

TEN PACK  
95362 | CPT: 29540

HALLUX CONTROL  
STRAP™



\*shown with Midfoot  
Compression Wrap  
(sold separately)

TEN PACK  
9534\*T | CPT: 29550

NIGHT SPLINT



3320\* | HCPCS: L4397

NON-PNEUMATIC  
WALKING BOOT



SHORT

TALL

SHORT  
3311\* | HCPCS: L4387  
TALL  
3312\* | HCPCS: L4387

PNEUMATIC WALKING  
BOOT



SHORT

TALL

SHORT  
3303\* | HCPCS: L4361  
TALL  
3310\* | HCPCS: L4361



HCPCS: A4467  
WITH COMPRESSION WRAP  
5351\* | HCPCS: A4467

TRILOK™



5351\* | HCPCS: L1902

VISCO ANKLE SKIN™



HCPCS: A4467  
WITH COMPRESSION WRAP  
5341\* | HCPCS: A4467

WEIL OSTEOTOMY  
STRAP™



\*shown with Midfoot  
Compression Wrap  
(sold separately)  
95372 | CPT: 29550

ANKLE SKIN™



KNEE

CROSSFIRE™ ALUMINUM



**PULL-ON**  
6186\* | HCPCS: L1820  
**WRAPAROUND**  
5186\* | HCPCS: L1820

CROSSFIRE™ OA



**WRAPAROUND**  
5190\* L or R | HCPCS: L1851

CROSSFIRE™ ROM



**PULL-ON**  
6180\* | HCPCS: L1833  
**PULL-ON WITH GEL**  
6187\* | HCPCS: L1833  
**WRAPAROUND**  
5180\* | HCPCS: L1833

CROSSFIRE™ CONFORMA



**PULL-ON**  
6185\* | HCPCS: L1820  
**PULL-ON WITH GEL**  
6188\* | HCPCS: L1820  
**WRAPAROUND**  
5185\* | HCPCS: L1820

GLADIATOR™



**PULL-ON (ULTIMA)**  
4134\* | HCPCS: L1833  
**PULL-ON (STRATUS)**  
6136\* | HCPCS: L1833  
**WRAPAROUND**  
5136\* | HCPCS: L1833

HINGED KNEE SKIN™



**PULL-ON (ULTIMA)**  
4130\* | HCPCS: L1820  
**PULL-ON (STRATUS)**  
6130\* | HCPCS: L1820  
**WRAPAROUND**  
5133\* | HCPCS: L1820

PATELLA STABILIZER W/CONFORMA™ HINGE



**ULTIMA**  
4132\* | HCPCS: L1820  
**STRATUS**  
6132\* | HCPCS: L1820  
**WITHOUT A HINGE**  
4125\* | HCPCS: A4467

PREMIUM J™ W/CONFORMA™ HINGE



**W/CONFORMA HINGE**  
6138\* L or R | HCPCS: L1820  
**WITHOUT A HINGE**  
6163\* L or R | HCPCS: A4467

Q BABY™



5165\* | HCPCS: A4467

Q BRACE™ W/CONFORMA™ HINGE



**W/CONFORMA HINGE PULL-ON**  
4136\* | HCPCS: L1820  
**WITHOUT A HINGE PULL-ON**  
4161\* | HCPCS: A4467  
**W/CONFORMA HINGE WRAPAROUND**  
5138\* | HCPCS: L1820  
**WITHOUT A HINGE WRAPAROUND**  
5163\* | HCPCS: A4467

Q LOK™ ROM



**PULL-ON**  
6166\* | HCPCS: L1833  
**WRAPAROUND**  
5166\* | HCPCS: L1833

KNEE SKIN™



4111\*O or C | HCPCS: A4467

TEN-7™



5146\* | HCPCS: A4467

UNDER SKIN™



4140\* | HCPCS: A4467

VISCO KNEE SKIN™ W/CONFORMA™ HINGE



**W/CONFORMA HINGE**  
6137\* | HCPCS: L1820  
**WITHOUT A HINGE**  
6171\* | HCPCS: A4467

SPINE

BACK FLEX™



3730\* | HCPCS: L0642

BACK PLUS™



3740\* | HCPCS: L0648

BACK SKIN™



4710\* | HCPCS: A4467

**W/ FOAM LUMBAR PAD**  
4720\* | HCPCS: A4467

**W/ FLEXIBLE LUMBAR SUPPORT**  
4734\* | HCPCS: A4467

**W/ FOAM PAD & FLEXIBLE SUPPORT**  
4724\* | HCPCS: A4467

BAJA™ 627



3750\* | HCPCS: L0642

BAJA™ 631



3770\* | HCPCS: L0648

BAJA™ 637



3760\* | HCPCS: L0650

VECTOR™



3720\* | HCPCS: A4467

COMPRESSION



**CALF SKIN™**  
4630\* | HCPCS: A4467



**COMPRESSION SHORTS**  
4620\* | HCPCS: A4467  
**WITH GROIN WRAP**  
4621\* | HCPCS: A4467



**THIGH SKIN™**  
4610\* | HCPCS: A4467  
**WITH CINCH STRAP**  
4611\* | HCPCS: A4467