

TriLok® Unequaled Versatility

The TriLok has unparalleled versatility, delivering targeted unloading to medial and lateral structures. Positioned medially, the FootLok™ Strap supports the plantar fascia or the posterior tibial tendon. Positioned laterally, the FootLok Strap supports and protects the lateral ligaments and the peroneal tendons. The TriLok is effective with almost any type of footwear or with no footwear at all.

Medial Unloading:

Plantar Fasciitis – Posterior Tibial Tendon Dysfunction – Deltoid Sprains – Subtalar Malalignment



Undersleeve

1 Compression & Proprioception



FootLok Strap

2 Medial Unloading & Stabilization



Stirrup Strap

3 Ankle & Subtalar Control

Lateral Unloading:

Lateral Sprains – Peroneal Tendon Damage – Ankle Arthritis



Undersleeve

1 Compression & Proprioception



FootLok Strap

2 Lateral Unloading & Stabilization



Stirrup Strap

3 Ankle & Subtalar Control

Sizing

		By shoe size				
		TriLok™				
Size	*PN	Men's	Women's	Europe	UK	
S	1	5 – 6.5	6 – 7.5	38 – 39	4 – 5.5	
M	2	7 – 9.5	8 – 11	40 – 42	6 – 8.5	
L	3	10 – 12.5		43 – 45	9 – 11.5	
XL	4	13 – 16		46 – 49	12 – 16	

		HCS - Hallux Control Strap				
Size	*PN	Men's	Women's	Europe	UK	
S	1	5 – 10	6 – 11	38 – 43	4 – 9	
L	3	10.5 – 16	11.5+	44 – 50	9.5 – 16	

		Circumference from heel to top of the foot – back to heel		
		AFTR Products		
Size	*PN	Inches	Centimeters	
XS-S	1	11 – 13	27.9 – 33	
M-L	3	12.5 – 14	31.8 – 35.6	
XL-XXL	5	13.5 – 15	34.3 – 38.1	

Products

TriLok	Product #
	5360*

AFTR DC (Dorsal Closure)	Product #
w/ Rigid Stirrup Straps	5386*
w/ Rigid Stirrup Straps & figure 8 Compression Wrap	5386*CW

AFTR with Gel	Product #
w/ Rigid Stirrup Straps	5383*
w/ Rigid Stirrup Straps & figure 8 Compression Wrap	5383*CW

AFTR (No gel pads)	Product #
w/ Rigid Stirrup Straps	5380*
w/ Rigid Stirrup Straps & figure 8 Compression Wrap	5383*CW

HCS (Hallux Control Strap)		
Size	Qty	Product #
1 (S)	10-Pack	95341T
3 (L)	10-Pack	95343T

Weil Osteotomy Strap		
Size	Qty	Product #
Universal	10-Pack	95372

Calibrated Midfoot Compression Wrap		
Size	Qty	Product #
Universal	10-Pack	95362

Note: To indicate size, replace * with corresponding number from the sizing chart.



240 East Hersey St., Suite 2 Ashland, OR 97520-5202
 t: 800-541-2455 f: 541-482-2119 CropperMed.com



Patent #: 7,115,105
 ©2015 Cropper Medical, Inc.
 03/2015 PN: 12352E



The only brace to offer
Medial & Lateral Unloading

TriLok[®]

External Ligament System™

Versatility:

Can be applied for medial or lateral pathologies

Control:

Two independent straps stabilize the midfoot and hindfoot

Comfort:

Low profile and allows normal ankle flexion and extension



INDICATIONS

- Posterior Tibial Tendon Dysfunction
- Plantar Fasciitis
- Peroneal Tendinopathy
- Ankle Sprains
- Subtalar and Ankle Arthritis

AFTR™ Ankle & Foot Trauma Rehabilitation

The AFTR braces combine stabilization with high-level compression unique to Bio Skin® materials. Following trauma, AFTR braces help:

- Reduce Pain
- Improve Proprioception
- Increase ROM
- Reduce Recovery Time

AFTR DC™ (Dorsal Closure)

- Accommodates changes in foot volume
- Easy to apply during early post-op phase
- Includes medial and lateral silicone gel pads behind the malleoli
- Trimmable in length

AFTR / AFTR with Gel™

- Simple Application
- Trimmable in length
- Ideal for patients with less swelling
- AFTR with Gel includes medial and lateral silicone gel pads behind the malleoli



INDICATIONS

- Peroneal / Post-Tib Tendon Repair
- Lateral Ankle Repair
- Fusions
 - Midfoot
 - Hindfoot
 - Ankle
- Total Ankle Replacement
- Ankle Sprains
- Bunionectomy
- Forefoot Procedures
- Hindfoot Arthritis

WEIL Digit Strapping System™

Soft-tissue balancing during critical, early post-op healing.

HCS™ (Hallux Control Strap™)



- **Splinting** - Maintains the hallux in a neutral position during critical soft-tissue healing after bunionectomy.
- **Stretching** - Stretches extensor tendons after primary healing to improve ROM and allow normal gait.
- **Strengthening** - Provides resistance for patients to strengthen FHL for strong toe-off.

Shown with Calibrated Midfoot Compression Wrap. Also compatible with AFTR, AFTR DC and AFTR with Gel. Patent pending.

WEIL Osteotomy Strap™



- Prevents dorsal drift and floating toes by holding the toe in a neutral or plantar-flexed position.
- May be worn inside a walking boot, post-op shoe, or normal footwear.
- May be used non-operatively to treat metatarsalgia or plantar plate tears.

Shown with Calibrated Midfoot Compression Wrap. Also compatible with AFTR, AFTR DC and AFTR with Gel. Patent pending.

Compression Wrap

CMC Wrap™ (Calibrated Midfoot Compression)



- Provides an attachment point for the Weil Digit Strapping System without covering the ankle.
- Made of hypo-allergenic, non-slip material to prevent migration. Easily fits inside normal footwear.