

BAJA 637

WITH RIGID ANTERIOR PANEL

L0637/L0650

The Baja 637 uses patented technology to deliver targeted compression for immediate relief from a wide variety of conditions causing low back pain.

- Single, trim-to-fit size
- Bilateral pulley system focuses compression on the tissues supporting the lumbar and thoracic regions of the spine
- Circumferential compression stabilizes the spine and limits painful movement in the sagittal plane
- Lateral panels add stability and limit movement in the coronal plane
- Belt provides circumferential compression to increase intracavitary pressure and relieve pressure on spinal discs
- Pulleys make it easy to apply firm compression even for people with weak hands or limited dexterity
- Lightweight, breathable fabrics keep users cool and comfortable
- Rigid anterior panel increases surface area for greater comfort and extra support.
- Includes two hot/cold gel pads that slide into the lumbar pad
- Heat moldable posterior panel for a customized fit



Posterior panel extends from coccyx to thoracic spine.



Rigid anterior panel maximizes support and compression.



Rigid lateral panels for coronal control.

PRODUCT

PART #

Baja 637 with Rigid Anterior Panel

3776**RAP



Pro.BioSkin.com

CALL: 800-541-2455 | FAX: 541-482-2119 | EMAIL: orders@bioskin.com